

# DOSTA LIČNO STOJNE

(Makedonia)

MUSIC: Cassette Tape AK-015, side A band 1

FORMATION: W hold, mixed line

---

METER: 7/8 (SQQ) 1 2 3

PATTERN

---

Meas

- Part 1 (singing part)
- 1-3 Lesnoto step on R fwd(ct.1); čukče on R (ct.2); Step on L fwd (ct.3) step on R to R, face ctr (ct 4), čukče on R twice, lift L knee up in front (ct.5,6) Repeat cts( 4-6) with opp ft. (cts.7-9)
- 4-24 Repeat meas 1-3, 7 more times

- Part 2
- 1-2 Repeat Part 1, Meas. 1-2
- 3 Step on L diag. L back (ct.1); Step on R to R and face LOD(cts.2,3)
- 4 Step on L in front of R (ct.1); Step on ~~L~~ R fwd(cts.2,3)
- 5-16 Repeat Meas. 1-4, 3 more times. *R+L*

## DOSTA LIČNO STOJNE

Dosta Dosta lično bre Stojne  
Sluginče sluguvaš  
//Sluginče sluguvaš bre Stojne  
Sekoj sboru zboruva//

Ajde ajde lično bre Stojne  
Begaj ti po mene  
//Begaj ti po mene bre Stojne  
Ja stradam za tebe//

Ne mozam ne mozam bre ludo  
Ništo jas ci nemam  
//Sirače Sum jas pora snalo  
Bez majka bez tatko//

Ništo ne ti sakam bre Stojne  
Za Sakajme od srce  
//Ljubov taija sakam bre Stojne  
Naj sakana do groba//